


How is GoCheck Kids Clinically Validated and How are Referring Thresholds Set?

GoCheck Kids uses eccentric photorefraction - Eccentric photorefraction uses the iPhone's flash that is eccentrically positioned relative to the camera aperture. A bright crescent appears in the pupillary reflex when the child's eye has sufficient hyperopia or myopia along the meridian of the flash eccentricity and the app performs estimation of refraction through proprietary algorithms and provides immediate viewing of the results.



Principles of Photorefraction



The photorefraction value is calculated using the delta crescent center methodology, where by the crescents are measured from the edge of the originating pupillary margin as they approach or cross the center of the pupil. The calculation is based on the position and size of the crescent, as well as the pupil size.

To clinically validate photoscreening the following takes place:

- 1.) Children who are being seen by pediatric ophthalmology are consented and the patient has an image taken with the app.
- 2.) That image is processed and a "photo refraction" value is determined.
- 3.) The child then receives a cycloplegic exam and the cycloplegic refraction is determined.
- 4.) The results of the cycloplegic exam are compared to the AAPOS criteria for referral (the pediatric ophthalmology group that has determined at which refractive level a child is most at risk for the development of amblyopia).

AAPOS Referral criteria:

Age, months	Astigmatism	Hyperopia	Anisometropia	Myopia
12-30	>2.0 D	>4.5 D	>2.5 D	>-3.5 D
31-48	>2.0 D	>4.0 D	>2.0 D	>-3.0 D
>48	>1.5 D	>3.5 D	>1.5 D	>-1.5 D

most children are slightly hyperopic, but to avoid confusion with primary care and non-ophthalmology providers, the report displays 00.0D. As shown:


**Cook Children's
Magnolia**



Visual Screening Report

Date of Report: December 3, 2018
ePhild: RS160603

Provider Name: Frank McGehee



No risk factors identified at this time

Test Date: 12/3/2018 Patient Age: 4
Taken with Device: FCGT64QTHG03

Refractive Data	Right	Left	Photorefraction Threshold*
Hyperopia (D)	0.00	0.00	1.75
Myopia (D)	0.00	0.00	2.00
Anisometropia (D)	0.00		1.25

Legend

XX.XX photorefraction in diopters (D), below risk threshold

XX.XX photorefraction in diopters (D), above risk threshold

*Photorefraction thresholds yield sensitivity and specificity compared to cycloplegic refraction using the 2013 AAPOS referral criteria


In a child with crescents present on the image, the crescents are measured using proprietary software and algorithms, the photorefraction is calculated and displayed as shown, in diopters of photorefraction:

Visual Screening Report

Date of Report: August 2, 2018

ePhid: MM180802

Provider Name:



Risk factors identified

Test Date: 8/2/2018 Patient Age: 4
Taken with Device: 2017-02-40

Refractive Data	Right	Left	Photorefraction Threshold*
Hyperopia (D)	1.87	1.38	1.24
Myopia (D)	0.00	0.00	2.00
Anisometropia (D)	0.49		0.70

Legend

XX.XX photorefraction in diopters (D), below risk threshold

XX.XX photorefraction in diopters (D), above risk threshold

*Photorefraction thresholds yield sensitivity and specificity compared to cycloplegic refraction using the 2013 AAPOS referral criteria.

The referring thresholds shown photorefraction thresholds that have been determined using the above methodology to yield the best sensitivity and specificity when calibrated using cycloplegic refraction.